

Technical Information

OUA260

Flow assembly for sensors OUSAFxx and OUSTF10



Application

The flow assembly OUA260 can be combined with a wide array of optical sensors. Depending on the sensor used, the flow assembly can be deployed in the following applications:

- Chromatography monitoring
- Filtration monitoring
- Color measurement
- Centrifuge control
- Measurement of protein concentrations
- Turbidity measurement

Your benefits


- Fast measured values ensure maximum product yield
- Low sample volume reduces product loss
- Flexible use owing to wide range of nominal diameters and process connections, such as Tri-Clamp, flanges, threads etc.
- Hygienic and rugged versions thanks to variety of materials for assembly, seals and windows
- Complies with the requirements of the Life Sciences industry:
Certified biocompatibility in relation to biological reactivity in accordance with USP <87> and USP <88> Cl. VI, FDA-listed seals and hygienic, electropolished surfaces (Ra = 0.38 µm (15 µinch))

Function and system design

Measuring system

A complete measuring system comprises:

- Liquiline CM44P transmitter
- Photometer sensor, e.g. OUSAF44
- Flow assembly OUA260
- Cable set CUK80



A0031510

1 Measuring system with OUA260

- 1 Pipe
- 2 Transmitter CM44P
- 3 Cable set CUK80
- 4 Sensor: detector
- 5 Flow assembly OUA260
- 6 Sensor: light source (lamp)
- 7 Cable set CUK80

Options

Precision optical path length (POPL) adjustment system for OUA260


This option for the flow assembly guarantees the precise setting of the distance between the windows. The path length adjustment system consists of adjustable window rings and a certified measuring gage that determines the precise distance between the windows. This option allows the operator to set the optical path length down to a minimum of length of 0.5 mm.

The advantages:

- Increased measuring range
- Better measured value reproducibility
- Consistent readings across different devices
- Direct concurrence with laboratory results

Combined with Easycal, it enables calibration traceable to NIST and eliminates the need for time-consuming calibration with liquid standards.

The optical windows can be purged with dry air or nitrogen through pneumatic ports, preventing the formation of condensate on the optical windows.



A0025475

2 Example of purge air supply

- 1 Compressed air or nitrogen supply
- 2 Air drier (not required for nitrogen)
- 3 Pressure regulator
- 4 Flow controller
- 5 Assembly OUA260

The purge gas must be clean and dry (ultra zero air).


Maximum pressure:	0.07 bar (1 psi)
Flow rate:	50 to 100 ml/min

Installation

Installation instructions

The assembly is available with a range of process connections. It can be installed either directly in a process line or in a bypass line.

- ▶ Make sure that the optical windows of the assembly are fully immersed in the medium.
- ▶ Avoid installation positions in which air bubbles can form.
- ▶ Install the flow assembly upstream from the pressure regulators.



A0028250

3 Mounting angles. The arrows indicate the direction of medium flow in the pipe.

- A Suitable installation position, better than C
- B Ideal, best installation position
- C Acceptable installation position
- D Installation position to be avoided
- E Unacceptable installation position

Process

Process temperature range and pressure range


The process temperature range and pressure range depend on the process connection, material and pipe size.


Process connection	Pipe size	Nominal pressure	Temperature
Tri-Clamp 316L	0.25 to 2"	16 bar (230 psi)	0 to 130 °C (32 to 266 °F)
Tri-Clamp 316L	2.5 to 4"	10 bar (150 psi)	0 to 130 °C (32 to 266 °F)
Tri-Clamp PVDF	0.25", 0.5", 0.75"	4 bar (58 psi)	0 to 130 °C (32 to 266 °F)
Flange ASME RF Class150, 316L	All	10 bar (150 psi)	0 to 130 °C (32 to 266 °F)
Flange RF Class150, 316L	All	20 bar (300 psi)	0 to 130 °C (32 to 266 °F)
NPT 316L	All	20 bar (300 psi)	0 to 130 °C (32 to 266 °F)
NPT PVDF, plastic fittings	All	4 bar (58 psi)	0 to 130 °C (32 to 266 °F)
NPT PVDF, metal fittings	All	2 bar (29 psi)	0 to 35 °C (32 to 95 °F)

 Please comply with the maximum permitted process temperature of the sensor.

Mechanical construction


Design



 4 Cross-sectional view

- 1 Process connection
- 2 Window ring
- 3 Connection thread for sensor housing (lamp)
- 4 Optical windows
- 5 Connection thread for sensor housing (detector)

Dimensions



A0024809

5 Dimensions of flow assembly OUA260


- A Flange spacing
 B internal diameter
 C Flange diameter

Process connection	Pipe size	A	B	C
Tri-Clamp	¼"	82.5 mm (3.25")	4.6mm (0.18")	25 mm (0.98")
Tri-Clamp	½"	82.5 mm (3.25")	9.4 mm (0.37")	25 mm (0.98")
Tri-Clamp	¾"	82.5 mm (3.25")	15.2 mm (0.60")	25 mm (0.98")
Tri-Clamp	1"	82.5 mm (3.25")	22.1 mm (0.87")	50.3 mm (1.98")
Tri-Clamp	1½"	82.5 mm (3.25")	36.1 mm (1.42")	50.3 mm (1.98")
Tri-Clamp	2"	82.5 mm (3.25")	47.2 mm (1.86")	64 mm (2.52")
Tri-Clamp	2½"	88.9 mm (3.50")	59.9 mm (2.36")	77.5 mm (3.05")
Tri-Clamp	3"	114.3 mm (4.50")	72.6 mm (2.86")	90.9 mm (3.58")
Tri-Clamp	4"	124.0 mm (4.88")	96.8 mm (3.81")	118.9 mm (4.68")
RFF150	1"	174.7 mm (6.88")	25.4 mm (1.00")	107.9 mm (4.25")
RFF150	2"	190.5 mm (7.50")	47.5 mm (1.87")	152.4 mm (6.00")
RFF150	3"	203.2 mm (8.00")	69.8 mm (2.75")	190.5 mm (7.50")
RFF150	4"	228.6 mm (9.00")	95.2 mm (3.75")	228.6 mm (9.00")
RFF300	1"	174.7 mm (6.88")	25.4 mm (1.00")	124.0 mm (4.88")
RFF300	2"	190.5 mm (7.50")	47.5 mm (1.87")	165.1 mm (6.50")
RFF300	3"	203.2 mm (8.00")	69.8 mm (2.75")	209.6 mm (8.25")
RFF300	4"	228.6 mm (9.00")	95.2 mm (3.75")	254.0 mm (10.00")
NPT-SS	½"	148.6 mm (5.85")	½" Standard NPT	N/A
NPT-SS	1"	101.6 mm (4.00")	1" Standard NPT	N/A
NPT-SS	2"	101.6 mm (4.00")	2" Standard NPT	N/A
NPT-PVDF	½"	71.1 mm (2.80")	½" Standard NPT	N/A
NPT-PVDF	1"	101.6 mm (4.00")	1" Standard NPT	N/A

Please contact the manufacturer for the sizes of Swagelock BVCO, Swagelock tube end and tube adapter.

Window types and path lengths

For both window types, the length is measured over the entire length.



A0024807

6 Length measurement of both window types

Example:

To obtain a path length of 10 mm with the Tri-Clamp 2.5" process connection, take one window with a length of 34 mm and one with a length of 36.8 mm.

Window types and path lengths for the various pipe sizes with the Tri-Clamp process connection

Path length	0.25" 0.50" 0.75"	1.0" LV 1.5" LV	2.0"	2.5"	3.0"	4.0"
0.5 mm POPL	19 + 18.5	24 + 23.5	33.5 + 34			
1 mm POPL	18 + 19	23 + 24	33.5 + 33.5			
2 mm	18 + 18	23 + 23				
2 mm POPL	18 + 18	23 + 23				
5 mm	14 + 19	19 + 24	31.5 31.5			
5 mm POPL	16.5 + 16.5	21.5 + 21.5	31.5 + 31.5			
10 mm	14 + 14	19 + 19	29 + 29	34 + 36.8		
20 mm	9 + 9	14 + 14	24 + 24	29 + 31.5	34 + 34	
30 mm		9 + 9	19 + 19	21.5 + 29	29 + 29	
40 mm			14 + 14	19 + 21.5	24 + 24	36.8 + 36.8
50 mm			9 + 9	14 + 16.5	19 + 19	31.5 + 31.5
60 mm				9 + 9	14 + 14	24 + 29
70 mm					9 + 9	21.5 + 21.5
80 mm						16.5 + 16.5
90 mm						9 + 14

Dimensions of the window types given in mm (e.g. 19 mm + 18.5 mm)

Window types and path lengths for the various pipe sizes with NPT SS and RFF 150/300 process connections

Path length	NPT SS 0.5" / 1.0" / 2.0"	RFF 150/300 1.0" / 2.0"	RFF 150/300 3.0"	RFF 150/300 4.0"
0.5 mm POPL	33.5 + 34	33.5 + 34		
1 mm POPL	33.5 + 33.5	33.5 + 33.5		
2 mm				
2 mm POPL				
5 mm	31.5 + 31.5	31.5 + 31.5		
5 mm POPL	31.5 + 31.5	31.5 + 31.5		
10 mm	29 + 29	29 + 29		

Path length	NPT SS 0.5" / 1.0" / 2.0"	RFF 150/300 1.0" / 2.0"	RFF 150/300 3.0"	RFF 150/300 4.0"
20 mm	24 + 24	24 + 24	34 + 34	
30 mm	19 + 19	19 + 19	29 + 29	
40 mm	14 + 14	14 + 14	24 + 24	36.8 + 36.8
50 mm	9 + 9	9 + 9	14 + 24	31.5 + 31.5
60 mm			14 + 14	24 + 29
70 mm			9 + 9	21.5 + 21.5
80 mm				16.5 + 16.5
90 mm				9 + 14

Dimensions of the window types given in mm (e.g. 19 mm + 18.5 mm)

Materials	Flow assembly:	Stainless steel 1.4404/1.4435 (AISI 316L), PVDF, other materials available on request
	Window:	Borosilicate, quartz, sapphire
	O-rings:	VITON-FDA, silicone FDA, EPDM-FDA, KALREZ-FDA
	PVDF is not suitable for all hazardous areas.	

Weight	Tri-Clamp ¼"	SS 316L:	1.14 kg (2.51 lbs)
	Tri-Clamp 1"	SS 316L:	1.39 kg (3.07 lbs)
	Tri-Clamp 2"	SS 316L:	1.88 kg (4.15 lbs)
	Tri-Clamp 4"	SS 316L:	3.38 kg (7.45 lbs)


Ordering information

Product page www.endress.com/oua260

Product Configurator

On the product page there is a "Configure" button to the right of the product image **Configure**.

1. Click this button.
 - ↳ The Configurator opens in a separate window.
2. Select all the options to configure the device in line with your requirements.
 - ↳ In this way, you receive a valid and complete order code for the device.
3. Export the order code as a PDF or Excel file. To do so, click the appropriate button on the right above the selection window.

 For many products you also have the option of downloading CAD or 2D drawings of the selected product version. Click the tab for this **CAD** and select the desired file type using picklists.

Scope of delivery

- The delivery comprises:
- Assembly in the version ordered
 - Operating Instructions

Available standard versions

The standard versions are listed in the following tables.

Process connections with nominal diameters available as standard:

Process connection (OUA260- *xx*****)		Nominal diameter (OUA260- ***x*****)										
		A 1/4"	B 3/8"	C 1/2"	D 3/4"	E 1" LV	F 1"	G 1½" LV	I 2"	J 2½"	K 3"	L 4"
A1	TRI-Clamp SS	✓		✓	✓	✓		✓	✓	✓	✓	✓
A2	Tri-Clamp PVDF	✓		✓	✓							
B1	Flange ASME RF Class 150						✓		✓		✓	✓
B2	Flange ASME RF Class 300						✓		✓		✓	✓
D1	Female NPT SS			✓			✓		✓			
D2	Female NPT PVDF			✓			✓					
E1	Swagelock BVCO	✓	✓	✓								
E2	Swagelock tube end	✓	✓	✓	✓	✓						
F1	Tube adapter		✓									

Path lengths with nominal diameters available as standard:

Path length (OUA260- ****xx*****)		Nominal diameter (OUA260- ***x*****)													
		A 1/4"	B 3/8"	C 1/2" (1)	C 1/2" (2)	C 1/2" (3)	D 3/4"	D 3/4" (3)	E 1" LV	E 1" Std	G 1½" LV	I 2"	J 2½"	K 3"	L 4"
01	0.5 mm /POPL	✓	✓	✓	✓		✓		✓		✓				
03	1 mm/POPL	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓			
04	2 mm	✓	✓	✓	✓		✓		✓		✓				
05	2 mm/POPL	✓	✓	✓	✓		✓		✓		✓				
06	5 mm	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓			
07	5 mm/POPL	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓			
08	10 mm	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		

Path length (OUA260-***x*****)		Nominal diameter (OUA260-***x*****)													
		A 1/4"	B 3/8"	C 1/2" (1)	C 1/2" (2)	C 1/2" (3)	D 3/4"	D 3/4" (3)	E 1" LV	E 1" Std	G 1 1/2" LV	I 2"	J 2 1/2"	K 3"	L 4"
09	20 mm	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
10	30 mm				✓	✓			✓	✓	✓	✓	✓	✓	✓
11	40 mm					✓			✓		✓	✓	✓	✓	✓
12	50 mm					✓			✓		✓	✓	✓	✓	✓
13	60 mm											✓	✓	✓	✓
14	70 mm													✓	✓
15	80 mm														✓
16	90 mm														✓

C 1/2" (1) = Tri-Clamp and Swagelok process connections (OUA260-***x*****) ; A1, A2, E1, E2)

C 1/2" (2) = PVDF-NPT process connection (OUA260-***x*****) ; D2)

C 1/2" (3) and D 3/4" (3) = SS-NPT process connection (OUA260-***x*****) ; D1)

Accessories

The following are the most important accessories available at the time this documentation was issued.

- ▶ For accessories not listed here, please contact your Service or Sales Center.

Sensors

OUSAF44

- Optical sensor for measuring UV absorption
- Variety of materials and process connections available
- Hygienic design
- Product Configurator on the product page: www.endress.com/ousaf44



Technical Information TI00416C

OUSAF12

- Optical sensor for the measurement of absorbance
- Variety of materials and process connections available
- Product Configurator on the product page: www.endress.com/ousaf12



Technical Information TI00497C

OUSAF22

- Optical sensor for measuring color concentrations
- Variety of materials and process connections available
- Product Configurator on the product page: www.endress.com/ousaf22



Technical Information TI00472C

OUSTF10

- Optical sensor for measuring turbidity and undissolved solids
- Variety of materials and process connections available
- Product Configurator on the product page: www.endress.com/oustf10



Technical Information TI00500C

OUSAF46

- Optical sensor for measuring UV absorption
- Two individually configurable measuring channels
- Product Configurator on the product page: www.endress.com/ousaf46



Technical Information TI01190C

www.addresses.endress.com
